



Tableau des valeurs de couple des instruments MICRO-MEGA®

Torque values for MICRO-MEGA® instruments

Drehmomentwerte für die MICRO-MEGA® Instrumente

MEM		N°	Conicité / Taper / Konizität	Zone de travail / Working area / Arbeitsgebiet	N.Cm	Vitesse (tr/min) / Speed (rpm) / Geschwindigkeit (U/min)
1	ENDOFLARE®	25	12%	1/4 LT/WL/AL	3	400
2	G-Files™ G1	12	3%	LT/WL/AL	1.2	400
	G-Files™ G2	17	3%	LT/WL/AL	1.2	400
3	One Shape®	25	6%	LT/WL/AL	2.5	400
4	One Shape® Apical 1	30	6%	LT/WL/AL	1	400
	One Shape® Apical 2	37	6%	LT/WL/AL	1	400
5	Revo-S™ SC1	25	6%	2/3 LT/WL/AL	0.8	400
	Revo-S™ SC2	25	4%	LT/WL/AL	0.8	400
	Revo-S™ SU	25	6%	LT/WL/AL	0.8	400
	Revo-S™ AS30	30	6%	LT/WL/AL	0.8	400
	Revo-S™ AS35	35	6%	LT/WL/AL	0.8	400
	Revo-S™ AS40	40	6%	LT/WL/AL	0.8	400

MM•control

